

# Summer project presentation

Closed loop neurofeedback prototype for adaptive brain stimulation


Student:

Zeca Buclet

Supervisor:


Professor Marie-Hélène Boudrias

# System overview



# Summer project overview

Objective: Extract MRBD and trigger stimulation if MRBD below threshold



Current Opinion in Biomedical Engineering

## Tasks

1. Feasibility study
2. Online feature extraction
3. ~~Trigger stimulation~~

# Feasibility Study



## Literature Review

- Brain stimulation and why use it in stroke recovery
- Task of interest & detectable features
- Neurofeedback and applications



## Requirements

- Real time data access
- Low latency transmission
- Easy integration and system translation

# Unicorn Recorder



Acquisition

Recording

Network

Device

## Trigger UDP Input Settings

Enabled ☒

IP

Port

## Raw Data UDP Output Settings

Enabled ☒

IP

Port

## Raw Data LSL Output Settings

Enabled ☐

Streamname

## Processed Data UDP Output Settings

Enabled ☒

IP

Port

## Processed Data LSL Output Settings


Enabled ☒

Streamname

# Feasibility study

Find appropriate data access protocol:

Option	Comments
BCPy2000	Outdated, complicated
NeurofeedbackLab	Matlab based, EEGLAB plugins
BrainVision RecView	Live processing, no control over data
LabStreamingLayer (LSL)	Modern, cross-compatible, Python API, low latency

 Chosen protocol: LabStreamingLayer

# Feasibility study: Tools used

## Hardware:

- Personal computer
- Unicorn Hybrid Black

### ✓ Pros

- Wireless
- Dry electrodes
- Easy to use

### ✗ Cons

- Noisier
- Lower sampling rate

## Software:

- Unicorn Suite Hybrid Black
- Python
- LabStreamingLayer




# Feature extraction: Record and log data

## Goal

Detect alpha power shifts in real time → Compare eyes open and eyes closed recording

## LSL protocol

- Data streamed from UnicornRecorderLSLStream (Processed data)
- Python script to access and log data into csv file










Image source: Unicorn Suite Hybrid Black – Unicorn Recorder Screenshot

# Feature extraction: Offline analysis alpha power

Alpha power change between tasks (eyes open/closed) relative to baseline



## ✗ Problems


- Data noisier than expected
- Big variability across trials
- Small baseline values  $\Rightarrow$  big changes

## ✓ Detection

- Clear alpha band increase

# Feature extraction: Online extraction

## 3) Alpha power extraction online logic:




### Pseudo code

```
LOOP while streaming is active:
    Pull_sample_from_stream()
    if enough data to compute on a new window:
        if bad_window:
            BREAK
        Compute_Welch_PSD()
        Compute_alpha_power()
        Log data to csv
END LOOP
```

Parameters	Value
Sample rate	250Hz
Window size	2s
Buffer size	Window x sample rate
Step size	10% buffer (90% overlap)

# Feature extraction: Online analysis alpha power



## ✗ Limitations

- Very noisy
- Large standard deviation

## ✓ Results

- Successful real-time detection of state changes (20ms)
- Validated pipeline for recording, plotting, and feature extraction

# Feature extraction: Offline analysis MRBD

## Goal


Detect movement-related beta desynchronization (MRBD)

## Method


- Extend alpha power pipeline to beta band
- 4 second gripping task — 6s rest  $\Rightarrow$  50 trials
- Morlet wavelet time frequency analysis over C3

# Feature extraction: Offline analysis MRBD

✓ MRBD Validation on provided data




- Clear 50% desynchronization.
- Post movement rebound




# Feature extraction: Offline analysis MRBD

✗ MRBD extraction my own data





- No clear desynchronization.
- Possible noise/artifact issues



# Feature extraction: Offline analysis MRBD

Movement vs No Movement comparison



## Interpretation

- Task related signal is inexistent
- Beta power baseline is too small and data is too **noisy** which creates very large ERDs





# Recap & Next steps

## ✓ Success

- Real-time feature extraction feasible (alpha)
- Mean alpha change computed live per task

## ⚠ Limitations


- MRBD detection harder than expected
- Higher-quality recordings required
- More robust analysis & protocol

## 🚀 Next steps

- Use BrainVision Cap for better data quality
- Optimize preprocessing (filtering, resampling)
- Improve experimental protocol




# Optimal EEG setup for closed loop system



Current Opinion in Biomedical Engineering

# Optimal EEG setup for closed loop system



Current Opinion in Biomedical Engineering

# Optimal EEG setup for closed loop system

## To Buy

- **Soterix MxN Pro** – fully software-controlled stimulation system

**Current limitation:** 4×1 HD-tDCS and 1×1 tES cannot be modulated via software.



1

## To Develop

- Live MRBD computation
- Stimulation trigger



2

**Images sources:** 1) <https://soterixmedical.com/static/images/mxn33/mxn-pro-main-3.png>

2) <https://upload.wikimedia.org/wikipedia/commons/c/c3/Python-logo-notext.svg>

# What I learned



## Technical skills

- Real time processing pipeline using Python and LSL
- Online and offline feature extraction
- Hands on experience with EEG hardware



## Research experience

- Learned importance of iterative testing and refining methods
- Problem solving on real use cases
- Patience and resilience
- Unexpected results are as important as expected ones



Thank you for your attention!

# References

## Bibliography:

I. Iturrate, M. Pereira, and J. del R. Millán, “Closed-loop electrical neurostimulation: Challenges and opportunities,” *Current Opinion in Biomedical Engineering*, vol. 8, pp. 28–37, 2018. doi: 10.1016/j.cobme.2018.09.007

Morales Fajardo K, Yan X, Lungoci G, Casado Sánchez M, Mitsis GD, Boudrias MH. The Modulatory Effects of Transcranial Alternating Current Stimulation on Brain Oscillatory Patterns in the Beta Band in Healthy Older Adults. *Brain Sci.* 2024 Dec 20;14(12):1284. doi: 10.3390/brainsci14121284

## Images:

- 1) <https://www.gtec.at/wp-content/uploads/2023/09/unicorn-hybrid-black-bundle.jpg>
- 2) <https://soterixmedical.com/static/images/mxn33/mxn-pro-main-3.png>
- 3) <https://upload.wikimedia.org/wikipedia/commons/c/c3/Python-logo-notext.svg>